





Aquatic areas accessible to everyone

The Hippocampe Pool makes it easier for **people with** reduced mobility to access aquatic and pool spaces.

Usable in water, easy to use and hygienic it accompanies users from changing rooms to pool **areas** while avoiding multiple transfers.



The Hippocampe allows access to your pool and aquatic complex







Accessibility

The Hippocampe allows people with reduced mobility to move around freely between changing rooms, showers, and pool areas.

It can be submerged in water to simplify access to swimming.



Hygiene and security

The seat and footrest cover is made of comfortable hydrophobic foam that protects its user.

easy maintenance, With Hippocampe can be rinsed and disinfected after each use.

FDA REGISTERED

included

Vipamat's goal is to insure a level of quality and performance that satisfies the requirements of the FDA.





Seat length - feet

28.5" - 72 cm



Back seat height

14.5" - 37 cm



Seat width

15.5"- 40 cm



Seat depth 14.5" - 37 cm

Total length

54" - 137 cm



Seat height - ground 13.5"- 34 cm



Interior width 19.5"- 50 cm



Total height 38" - 96 cm



Total width 27" - 68 cm



Weight

Weight: 11.5 kg - 25,3 lb Max user weight: 130 kg / 286 lb



Technical data

Seat : Comfortable hydrophobic foam (3cm thickness) covered by black fabric which dries quickly and avoids skin problems.

Frame: Stainless steel and anodized aluminum with a foam cover that is comfortale to the touch even in case of heat and cold exposure.

Push bar: Anodized aluminium tube with a handle covered in foam that improves comfort for the companion.

FDA REGISTERED

Vipamat's goal is to insure a level of quality and performance that satisfies the requirements of the FDA.





Join us!







www.vipamat.com

VIPAMAT
3 Rue Gustave Eiffel
56270 PLŒMEUR - FRANCE
Tél.: +33 2 97 86 24 87
Fax : +33 2 97 86 73 81
www.vipamat.com

contact@vipamat.com

USA/CANADA

VIPAMAT 3302 Canal St. Suite 62 Houston, TX 77003 Office: 832 516 8125 Mobile: 865 385 4350 contact-usa@vipamat.com